



CHEMGENE HLD₄I

Instrument Disinfectant Wipes

CHEMGENE HLD₄I disinfectant wipes are manufactured with the latest innovation in nonwoven material, offering unrivalled strength, absorbency and low lint properties at a very cost effective price. The wipes are impregnated with a multi-active disinfectant solution that is alcohol free and tested to EN standards as effective against a range of bacteria, viruses, fungi, yeasts and mycobacteria.

BENEFITS OF USE

ALCOHOL FREE

Compliant


- CHEMGENE HLD₄I wipes are manufactured in accordance with Quality Management System accredited to ISO 13485.
- Controlled dosage ensures effective disinfection with each wipe
- Wipes are tested to a range of European Standards
- Mechanical action enhances disinfection.

Convenient

- No requirement for pre-mixing, reducing wastage and risk of errors.
- Canister can be stored on a work top or in a wall mounted bracket.

Cost Effective

- Saves valuable time and reduces the need for cloths and cleaning solutions,
- 100 wipes per canister means low cost per wipe.
- Fully recyclable material and packaging



Class IIb registered

AREAS OF USE

These wipes are designed for use on high risk areas and on most surfaces, including glass, metal ceramics, plastic and rubber. The tub can be stored in the area of use for easy access to quickly disinfect surfaces at risk of contamination.

TECHNICAL INFORMATION

CHEMGENE HLD₄I disinfectant wipes are developed and tested by qualified chemists and microbiologists. They are then subjected to stringent EN testing procedures to ensure efficacy against a wide range of bacteria, viruses and fungi. Certificates are available on request.

PURPOSE

For the cleaning and disinfecting of non-surgically invasive devices, such as transvaginal and transrectal probes/transducers and for hard surfaces of non-invasive, non-lumened medical devices such as monitors, radiographic and diagnostic equipment.

Not to be used to disinfect surgically invasive medical devices.

MICROBIOLOGICAL EFFICACY

CHEMGENE HLD₄I has been tested effective against:

FUNGICIDAL ACTIVITY <i>Trichophyton rubrum</i> <i>Candida albicans</i>	BACTERICIDAL ACTIVITY <i>Pseudomonas aeruginosa</i> <i>Escherichia coli</i> <i>Staphylococcus aureus</i> <i>Enterococcus hirae</i> <i>Legionella pneumophila</i> <i>Salmonella typhimurium</i> <i>Listeria monocytogenes</i> Methicillin-resistant <i>Staphylococcus aureus</i> Methicillin-resistant <i>Staphylococcus pseudointermedius</i>
VIRUCIDAL ACTIVITY Hepatitis B Herpes Simplex Norovirus Vaccinia Adenovirus	
BACTERIOPHAGE ACTIVITY <i>Lactococcus lacti</i>	



CODE	PACK SIZE	CASE QUANTITY
XTM353	100 wipes	10



USE BIOCIDES SAFELY

**ALWAYS READ THE LABEL
AND PRODUCT INFORMATION
BEFORE USE**



www.medimark.myshopify.com

MEDIMARK
S C I E N T I F I C

The details given in this leaflet are correct at the time of going to press. The company reserves the right to change and improve products without prior notice.

Manufactured in the UK

For more information on our products or for technical support contact:

East Point, Seal, Kent, TN15 0EG, UK

T 08452 223344 F 01732 769530

info@medi-mark.co.uk www.medi-mark.co.uk

TRUSTED INFECTION CONTROL